




SELV mouse mAb

Catalog No	YP-mAb-11087
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SELV
Immunogen	Synthesized peptide derived from human SELV AA range: 250-300
Specificity	This antibody detects endogenous levels of SELV at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SELV
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Testis specific.
Function	function:May be involved in a redox-related process .,similarity:Belongs to the SelWTH family.,tissue specificity:Testis specific.,
Background	This gene encodes a selenoprotein, which contains a selenocysteine (Sec) residue at its active site. The selenocysteine is encoded by the UGA codon that normally signals translation termination. The 3' UTR of selenoprotein genes have a common stem-loop structure, the sec insertion sequence (SECIS), that is necessary for the recognition of UGA as a Sec codon rather than as a stop signal. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using SELV mouse mAb