




GLYL1 mouse mAb

Catalog No	YP-mAb-11095
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GLYATL1 GNAT
Immunogen	Synthesized peptide derived from human GLYL1 AA range: 42-92
Specificity	This antibody detects endogenous levels of GLYL1 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GLYL1;GLYATL1 GNAT
Observed Band	
Calculated Molecular Weight	
Cell Pathway	mitochondrion,
Tissue Specificity	Expressed in liver and kidney and, at lower levels, in pancreas, testis, ovary and stomach.
Function	catalytic activity:Acyl-CoA + glycine = CoA + N-acylglycine., function:Mitochondrial acyltransferase which transfers the acyl group to the N-terminus of glycine. Can conjugate a multitude of substrates to form a variety of N-acylglycines.,similarity:Belongs to the glycine N-acyltransferase family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using GLYL1 mouse mAb