




## NPPL3 mouse mAb

|                             |                                                                                                                                   |
|-----------------------------|-----------------------------------------------------------------------------------------------------------------------------------|
| Catalog No                  | YP-mAb-11111                                                                                                                      |
| Isotype                     | IgG                                                                                                                               |
| Reactivity                  | Human;Mouse;Rat                                                                                                                   |
| Applications                | WB                                                                                                                                |
| Gene Name                   | NPIPL3 KIAA0220                                                                                                                   |
| Immunogen                   | Synthesized peptide derived from human NPPL3 AA range: 831-881                                                                    |
| Specificity                 | This antibody detects endogenous levels of NPPL3 at Human                                                                         |
| Formulation                 | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                           |
| Source                      | Monoclonal, Mouse,IgG                                                                                                             |
| Purification                | The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.              |
| Dilution                    | WB 1:500-1:2000                                                                                                                   |
| Concentration               | 1 mg/ml                                                                                                                           |
| Purity                      | ≥90%                                                                                                                              |
| Storage Stability           | -20°C/1 year                                                                                                                      |
| Synonyms                    | NPPL3;NPIPL3 KIAA0220                                                                                                             |
| Observed Band               |                                                                                                                                   |
| Calculated Molecular Weight |                                                                                                                                   |
| Cell Pathway                | Membrane ; Single-pass membrane protein .                                                                                         |
| Tissue Specificity          |                                                                                                                                   |
| Function                    | caution:Could be the product of a pseudogene.,similarity:Belongs to the NPIP family.,                                             |
| Background                  |                                                                                                                                   |
| matters needing attention   | Avoid repeated freezing and thawing!                                                                                              |
| Usage suggestions           | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel. |



## Products Images



Western Blot analysis of various cells using NPPL3 mouse mAb