




LRC57 mouse mAb

Catalog No	YP-mAb-11115
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LRRC57
Immunogen	Synthesized peptide derived from human LRC57 AA range: 87-137
Specificity	This antibody detects endogenous levels of LRC57 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LRC57;LRRC57
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Lipid-anchor .
Tissue Specificity	
Function	similarity:Contains 7 LRR (leucine-rich) repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using LRC57 mouse mAb