




TLCD1 mouse mAb

Catalog No	YP-mAb-11123
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TLCD1
Immunogen	Synthesized peptide derived from human TLCD1 AA range: 162-212
Specificity	This antibody detects endogenous levels of TLCD1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TLCD1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	similarity:Contains 1 TLC (TRAM/LAG1/CLN8) domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using TLCD1 mouse mAb