




RNAS9 mouse mAb

Catalog No	YP-mAb-11130
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RNASE9
Immunogen	Synthesized peptide derived from human RNAS9 AA range: 86-136
Specificity	This antibody detects endogenous levels of RNAS9 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RNAS9;RNASE9
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	At the mRNA level, widely expressed (PubMed:18992174). At protein level, restricted to epididymis (PubMed:19137000). Expressed in spermatozoa (sperm head and neck), with higher levels on ejaculated and epididymal sperm than on testicular sperm (at protein level). Expressed in the epithelial cells of the epididymal tubule (at protein level). Not detected in muscle.
Function	similarity:Belongs to the pancreatic ribonuclease family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using RNAS9 mouse mAb