




PEBP4 mouse mAb

Catalog No	YP-mAb-11137
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PEBP4 CORK1 UNQ1933/PRO4408
Immunogen	Synthesized peptide derived from human PEBP4 AA range: 132-182
Specificity	This antibody detects endogenous levels of PEBP4 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PEBP-4;PEBP4;PEBP4 CORK1 UNQ1933/PRO4408
Observed Band	25kD
Calculated Molecular Weight	25kD
Cell Pathway	Secreted .
Tissue Specificity	
Function	function:Seems to promote cellular resistance to TNF-induced apoptosis by inhibiting activation of the Raf-1/MEK/ERK pathway, JNK and phosphatidylethanolamine externalization.,similarity:Belongs to the phosphatidylethanolamine-binding protein family.,tissue specificity:Ubiquitously expressed. Highly expressed in tumor cells.,
Background	The phosphatidylethanolamine (PE)-binding proteins, including PEBP4, are an evolutionarily conserved family of proteins with pivotal biologic functions, such as lipid binding and inhibition of serine proteases (Wang et al., 2004 [PubMed 15302887]).[supplied by OMIM, Dec 2008],
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4° over night