




TPT2L mouse mAb

Catalog No	YP-mAb-11138
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TPTE2P1
Immunogen	Synthesized peptide derived from human TPT2L AA range: 40-90
Specificity	This antibody detects endogenous levels of TPT2L at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TPT2L;TPTE2P1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using TPT2L mouse mAb