




FA86C mouse mAb

Catalog No	YP-mAb-11140
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FAM86C1 FAM86C
Immunogen	Synthesized peptide derived from human FA86C AA range: 7-57
Specificity	This antibody detects endogenous levels of FA86C at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FA86C;FAM86C1 FAM86C
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	similarity:Belongs to the FAM86 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using FA86C mouse mAb