




HEAT6 mouse mAb

Catalog No	YP-mAb-11154
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	HEATR6 ABC1
Immunogen	Synthesized peptide derived from human HEAT6 AA range: 579-629
Specificity	This antibody detects endogenous levels of HEAT6 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HEAT6;HEATR6 ABC1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Amplified in breast cancer cell lines MCF-7 and BT-474.
Function	function:Amplification-dependent oncogene.,miscellaneous:17q23 region is one of the most commonly amplified regions in breast cancer and therefore may harbor genes important for breast cancer development and progression.,similarity:Contains 4 HEAT repeats.,tissue specificity:Amplified in breast cancer cell lines MCF7 and BT474.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using HEAT6 mouse mAb