




SPAC7 mouse mAb

Catalog No	YP-mAb-11175
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SPACA7 C13orf28
Immunogen	Synthesized peptide derived from human SPAC7 AA range: 115-165
Specificity	This antibody detects endogenous levels of SPAC7 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SPAC7;SPACA7 C13orf28
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted . Cytoplasmic vesicle, secretory vesicle, acrosome . Cytoplasmic vesicle, secretory vesicle, acrosome lumen . Localized in perinuclear pro-acrosomal granules in round spermatides. Localized between the inner and outer acrosomal membranes (matrix or lumen) in spermatozoa. Secreted during acrosome exocytosis. .
Tissue Specificity	Expressed in spermatozoa (at protein level) (PubMed:22495889).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using SPAC7 mouse mAb