



FA73A mouse mAb

Catalog No	YP-mAb-11177
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FAM73A
Immunogen	Synthesized peptide derived from human FA73A AA range: 36-86
Specificity	This antibody detects endogenous levels of FA73A at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FA73A;FAM73A
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Mitochondrion outer membrane ; Single-pass membrane protein .
Tissue Specificity	
Function	caution:It is uncertain whether Met-1 or Met-33 is the initiator.,similarity:Belongs to the FAM73 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

Western Blot analysis of various cells using FA73A mouse mAb

