




CG029 mouse mAb

Catalog No	YP-mAb-11179
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	C7orf29
Immunogen	Synthesized peptide derived from human CG029 AA range: 154-204
Specificity	This antibody detects endogenous levels of CG029 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CG029;C7orf29
Observed Band	
Calculated Molecular Weight	
Cell Pathway	nucleus,
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using CG029 mouse mAb