




FA21A mouse mAb

Catalog No	YP-mAb-11198
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FAM21A
Immunogen	Synthesized peptide derived from human FA21A AA range: 839-889
Specificity	This antibody detects endogenous levels of FA21A at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FA21A;FAM21A
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Early endosome membrane . Cell membrane .
Tissue Specificity	
Function	similarity:Belongs to the FAM21 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using FA21A mouse mAb