




CU024 mouse mAb

Catalog No	YP-mAb-11203
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LINC00114 C21orf24 NCRNA00114
Immunogen	Synthesized peptide derived from human CU024 AA range: 4-54
Specificity	This antibody detects endogenous levels of CU024 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CU024;LINC00114 C21orf24 NCRNA00114
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Expressed in a range of cell lines, including B-cell lymphoma and prostate.
Function	caution:Product of a dubious CDS prediction. Probable non-coding RNA. At least 6 alternative mRNA transcripts exist, with consensus splice donor and acceptor sites, but no coding potential nor murine ortholog. This protein sequence corresponds to the so-called isoforms C21orf24-6 and C21orf24-7., tissue specificity:Expressed in a range of cell lines, including B-cell lymphoma and prostate.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using CU024 mouse mAb