




PKHA3 mouse mAb

Catalog No	YP-mAb-11241
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PLEKHA3 FAPP1
Immunogen	Synthesized peptide derived from human PKHA3 AA range: 189-239
Specificity	This antibody detects endogenous levels of PKHA3 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PKHA3;PLEKHA3 FAPP1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Golgi apparatus, trans-Golgi network membrane ; Peripheral membrane protein . Localizes to ER-trans Golgi network (TGN) contacts sites. .
Tissue Specificity	Widely expressed.
Function	function: Binds specifically to phosphatidylinositol-4-phosphate (PtdIns4P), but not to other phosphoinositides.,similarity: Contains 1 PH domain.,tissue specificity: Widely expressed.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using PKHA3 mouse mAb