




SH3Y1 mouse mAb

Catalog No	YP-mAb-11244
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SH3YL1
Immunogen	Synthesized peptide derived from human SH3Y1 AA range: 10-60
Specificity	This antibody detects endogenous levels of SH3Y1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SH3Y1;SH3YL1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	ruffle membrane,
Tissue Specificity	
Function	similarity:Belongs to the SH3YL1 family.,similarity:Contains 1 SH3 domain., subunit:Interacts with SH3D19.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using SH3Y1 mouse mAb