




DDTL mouse mAb

Catalog No	YP-mAb-11248
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DDTL
Immunogen	Synthesized peptide derived from human DDTL AA range: 47-97
Specificity	This antibody detects endogenous levels of DDTL at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DDTL
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm .
Tissue Specificity	
Function	function:May have lyase activity .,similarity:Belongs to the MIF family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using DDTL mouse mAb