




ABCAB mouse mAb

Catalog No	YP-mAb-11254
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ABCA11P ABCA11
Immunogen	Synthesized peptide derived from human ABCAB AA range: 45-95
Specificity	This antibody detects endogenous levels of ABCAB at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ABCAB;ABCA11P ABCA11
Observed Band	
Calculated Molecular Weight	
Cell Pathway	integral component of membrane,
Tissue Specificity	
Function	caution:In contrast to other members of the family, it is much shorter and lacks ABC transporter domains. Defined as a pseudogene by HGNC. However, proteomic data on a specific peptide from PubMed:17924679 suggest the existence of this protein.,similarity:Belongs to the ABC transporter family. ABCA subfamily.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using ABCAB mouse mAb