



CG065 mouse mAb

Catalog No	YP-mAb-11258
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	C7orf65
Immunogen	Synthesized peptide derived from human CG065 AA range: 28-78
Specificity	This antibody detects endogenous levels of CG065 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CG065;C7orf65
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

Western Blot analysis of various cells using CG065 mouse mAb

