




SHSA3 mouse mAb

Catalog No	YP-mAb-11261
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SHISA3
Immunogen	Synthesized peptide derived from human SHSA3 AA range: 126-176
Specificity	This antibody detects endogenous levels of SHSA3 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SHSA3;SHISA3
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Endoplasmic reticulum membrane ; Single-pass type I membrane protein .
Tissue Specificity	
Function	function:Plays an essential role in the maturation of presomitic mesoderm cells by individual attenuation of both FGF and WNT signaling.,similarity:Belongs to the shisa protein family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using SHSA3 mouse mAb