




TM201 mouse mAb

Catalog No	YP-mAb-11262
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TMEM201 NET5 SAMP1
Immunogen	Synthesized peptide derived from human TM201 AA range: 290-340
Specificity	This antibody detects endogenous levels of TM201 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TM201;TMEM201 NET5 SAMP1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	[Isoform 2]: Nucleus inner membrane ; Multi-pass membrane protein . Cytoplasm, cytoskeleton, spindle pole . The C-terminal of isoform 2 is located on the nucleoplasmic side. During interphase, isoform 2 is distributed in the inner nuclear membrane in distinct micro-domains and during mitosis, it is found in the ER but it also localizes to the polar regions of the mitotic spindle. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using TM201 mouse mAb