




ZNF73 mouse mAb

Catalog No	YP-mAb-11270
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZNF73 ZNF186
Immunogen	Synthesized peptide derived from human ZNF73 AA range: 114-164
Specificity	This antibody detects endogenous levels of ZNF73 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZNF73;ZNF73 ZNF186
Observed Band	
Calculated Molecular Weight	
Cell Pathway	nucleus,
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using ZNF73 mouse mAb