




SMBT2 mouse mAb

Catalog No	YP-mAb-11275
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SFMBT2 KIAA1617
Immunogen	Synthesized peptide derived from human SMBT2 AA range: 156-206
Specificity	This antibody detects endogenous levels of SMBT2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SMBT2;SFMBT2 KIAA1617
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	
Function	similarity:Contains 1 SAM (sterile alpha motif) domain.,similarity:Contains 4 MBT repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using SMBT2 mouse mAb