




SCRN1 mouse mAb

Catalog No	YP-mAb-11291
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SCRN1 KIAA0193
Immunogen	Synthesized peptide derived from human SCR1 AA range: 164-214
Specificity	This antibody detects endogenous levels of SCR1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SES1;Secernin-1;secernin 1;KIAA0193;SCR1;SCR1 KIAA0193
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm .
Tissue Specificity	
Function	function:Regulates exocytosis in mast cells. Increases both the extent of secretion and the sensitivity of mast cells to stimulation with calcium., miscellaneous:'Secern' is an archaic English term meaning 'secrete'. , similarity:Belongs to the peptidase C69 family. Secernin subfamily.,
Background	This gene likely encodes a member of the secernin family of proteins. A similar protein in rat functions in regulation of exocytosis in mast cells. Alternatively spliced transcript variants have been described. [provided by RefSeq, Mar 2009],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using SCR11 mouse mAb