




LRRN3 mouse mAb

Catalog No	YP-mAb-11292
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LRRN3 Nbla10363 UNQ194/PRO220
Immunogen	Synthesized peptide derived from human LRRN3 AA range: 225-275
Specificity	This antibody detects endogenous levels of LRRN3 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LRRN3;LRRN3 Nbla10363 UNQ194/PRO220
Observed Band	
Calculated Molecular Weight	78kD
Cell Pathway	Membrane ; Single-pass type I membrane protein .
Tissue Specificity	
Function	similarity:Contains 1 fibronectin type-III domain.,similarity:Contains 1 Ig-like C2-type (immunoglobulin-like) domain.,similarity:Contains 12 LRR (leucine-rich) repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using LRRN3 mouse mAb