




ZDHC6 mouse mAb

Catalog No	YP-mAb-11294
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZDHHC6 ZNF376
Immunogen	Synthesized peptide derived from human ZDHC6 AA range: 160-210
Specificity	This antibody detects endogenous levels of ZDHC6 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1 :500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZDHC6;ZDHHC6 ZNF376
Observed Band	45kD
Calculated Molecular Weight	45kD
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein . When not palmitoylated, accumulates to dot-like structures in the endoplasmic reticulum (PubMed:28826475). .
Tissue Specificity	
Function	catalytic activity:Palmitoyl-CoA + protein-cysteine = S-palmitoyl protein + CoA., domain:The DHHC domain is required for palmitoyltransferase activity., similarity:Belongs to the DHHC palmitoyltransferase family.,similarity:Contains 1 DHHC-type zinc finger.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from HuvEc cells, primary antibody was diluted at 1:1000, 4° over night