




MOSC1 mouse mAb

Catalog No	YP-mAb-11302
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MARC1 MOSC1
Immunogen	Synthesized peptide derived from human MOSC1 AA range: 156-206
Specificity	This antibody detects endogenous levels of MOSC1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MARC1;MOSC1;MTARC1;MARC1 MOSC1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Mitochondrion outer membrane ; Single-pass type II membrane protein . Membrane ; Lipid-anchor . Mitochondrial import is mediated by AA 1-40 and requires ATP.
Tissue Specificity	
Function	cofactor:Molybdenum.,function:Probable oxidoreductase.,similarity:Contains 1 MOSC domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using MOSC1 mouse mAb