




LITD1 mouse mAb

Catalog No	YP-mAb-11307
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	L1TD1 ECAT11
Immunogen	Synthesized peptide derived from human LITD1 AA range: 797-847
Specificity	This antibody detects endogenous levels of LITD1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ECAT11;ES cell associated protein 11;ES cell-associated protein 11;LINE-1 type transposase domain-containing protein 1;LITD1;L1TD1 ECAT11
Observed Band	
Calculated Molecular Weight	
Cell Pathway	intracellular ribonucleoprotein complex,
Tissue Specificity	
Function	sequence caution:Translated as Leu.,sequence caution:Translation N-terminally extended.,similarity:Belongs to the transposase 22 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using LITD1 mouse mAb