




GAB4 mouse mAb

Catalog No	YP-mAb-11326
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GAB4
Immunogen	Synthesized peptide derived from human GAB4 AA range: 331-381
Specificity	This antibody detects endogenous levels of GAB4 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GRB2-associated binder 4;GRB2-associated binder 2-like;GRB2 associated binder 4;GRB2;GAB4
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	similarity:Belongs to the GAB family.,similarity:Contains 1 PH domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using GAB4 mouse mAb