




DPPA2 mouse mAb

Catalog No	YP-mAb-11353
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DPPA2 PESCRG1
Immunogen	Synthesized peptide derived from human DPPA2 AA range: 138-188
Specificity	This antibody detects endogenous levels of DPPA2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ECSA;Pluripotent embryonic stem cell-related gene 1 protein;PESCRG1; Developmental;DPPA2;DPPA2 PESCRG1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in embryonic stem cells. No expression is seen in 5 months embryo, mesenchymal stem cells, embryonic fibrocytes and adult tissues.
Function	function:May play a role in maintaining cell pluripotentiality.,similarity:Contains 1 SAP domain.,tissue specificity:Expressed in embryonic stem cells. No expression is seen in 5 months embryo, mesenchymal stem cells, embryonic fibrocytes and adult tissues.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using DPPA2 mouse mAb