




MEX3D mouse mAb

Catalog No	YP-mAb-11399
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MEX3D KIAA2031 RKHD1 RNF193
Immunogen	Synthesized peptide derived from human MEX3D AA range: 121-171
Specificity	This antibody detects endogenous levels of MEX3D at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MEX3D;MEX3D KIAA2031 RKHD1 RNF193
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Nucleus . Predominantly expressed in the cytoplasm and shuttles between the cytoplasm and the nucleus through the CRM1 export pathway. .
Tissue Specificity	Ubiquitously expressed in all the cell lines and tissues tested.
Function	domain:Binds RNA through its KH domains.,function:RNA binding protein, may be involved in post-transcriptional regulatory mechanisms.,PTM:Phosphorylated.,similarity:Contains 1 KH domain.,similarity:Contains 1 RING-type zinc finger.,similarity:Contains 2 KH domains.,subcellular location:Predominantly expressed in the cytoplasm and shuttles between the cytoplasm and the nucleus through the CRM1 export pathway.,tissue specificity:Ubiquitously expressed in all the cell lines and tissues tested.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using MEX3D mouse mAb