




DC4L2 mouse mAb

Catalog No	YP-mAb-11417
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DCAF4L2 WDR21C
Immunogen	Synthesized peptide derived from human DC4L2 AA range: 278-328
Specificity	This antibody detects endogenous levels of DC4L2 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DCAF4L2;WD repeat domain 21C;WDR21C;DC4L2;DCAF4L2 WDR21C
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using DC4L2 mouse mAb