




HECD2 mouse mAb

Catalog No	YP-mAb-11424
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	HECTD2
Immunogen	Synthesized peptide derived from human HECD2 AA range: 505-555
Specificity	This antibody detects endogenous levels of HECD2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HECT domain containing 2;HEC;HECD2;HECTD2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	nucleus,cytoplasm,cytosol,
Tissue Specificity	
Function	function:Probable E3 ubiquitin-protein ligase which accepts ubiquitin from an E2 ubiquitin-conjugating enzyme in the form of a thioester and then directly transfers the ubiquitin to targeted substrates.,pathway:Protein modification; protein ubiquitination.,similarity:Contains 1 HECT (E6AP-type E3 ubiquitin-protein ligase) domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using HECD2 mouse mAb