




KXDL1 mouse mAb


Catalog No	YP-mAb-11454
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA;IHC
Gene Name	KXD1 C19orf50
Immunogen	Synthesized peptide derived from human KXDL1 AA range: 34-84
Specificity	This antibody detects endogenous levels of KXDL1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000;IHC-p 1:50-300; ELISA 2000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KXD1;UPF0459 protein C19orf50;MSTP096;MST096;KxDL motif-containing protein 1;KXDL1;KXD1 C19orf50
Observed Band	19kD
Calculated Molecular Weight	19kD
Cell Pathway	Lysosome membrane .
Tissue Specificity	
Function	similarity:Belongs to the UPF0459 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from U2OS cells, primary antibody was diluted at 1:1000, 4° over night



Immunohistochemical analysis of paraffin-embedded human Gastric adenocarcinoma. 1. Antibody was diluted at 1:200(4° overnight). 2. Tris-EDTA, pH9.0 was used for antigen retrieval. 3. Secondary antibody was diluted at 1:200(room temperature, 45min).