



THAP6 mouse mAb

Catalog No	YP-mAb-11475
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	THAP6
Immunogen	Synthesized peptide derived from human THAP6 AA range: 39-89
Specificity	This antibody detects endogenous levels of THAP6 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	THAP6
Observed Band	
Calculated Molecular Weight	
Cell Pathway	microtubule cytoskeleton,
Tissue Specificity	
Function	similarity:Contains 1 THAP-type zinc finger.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

Western Blot analysis of various cells using THAP6 mouse mAb

