




TSN12 mouse mAb

Catalog No	YP-mAb-11490
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TSPAN12 NET2 TM4SF12 UNQ774/PRO1568
Immunogen	Synthesized peptide derived from human TSN12 AA range: 221-271
Specificity	This antibody detects endogenous levels of TSN12 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NET 2;NET2;Tetraspan NET 2;tetraspanin 12;TM4SF12;Tspan 12;TSPAN12;TSN12;TSPAN12 NET2 TM4SF12 UNQ774/PRO1568
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	similarity:Belongs to the tetraspanin (TM4SF) family.,
Background	The protein encoded by this gene is a member of the transmembrane 4 superfamily, also known as the tetraspanin family. Most of these members are cell-surface proteins that are characterized by the presence of four hydrophobic domains. The proteins mediate signal transduction events that play a role in the regulation of cell development, activation, growth and motility. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using TSN12 mouse mAb