




PL8L1 mouse mAb

Catalog No	YP-mAb-11510
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PLAC8L1
Immunogen	Synthesized peptide derived from human PL8L1 AA range: 40-90
Specificity	This antibody detects endogenous levels of PL8L1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PLAC8 like 1;PLAC8 like protein 1;PLAC8L1;PL8L1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	similarity:Belongs to the cornifelin family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using PL8L1 mouse mAb