




KBTB2 mouse mAb

Catalog No	YP-mAb-11516
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	KBTBD2 BKLHD1 KIAA1489 CGI-73
Immunogen	Synthesized peptide derived from human KBTB2 AA range: 359-409
Specificity	This antibody detects endogenous levels of KBTB2 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BKLHD1;KBTBD2;KIAA1489;KBTB2;KBTBD2 BKLHD1 KIAA1489 CGI-73
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cul3-RING ubiquitin ligase complex,
Tissue Specificity	
Function	similarity:Contains 1 BTB (POZ) domain.,similarity:Contains 5 Kelch repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using KBTB2 mouse mAb