




LRC59 mouse mAb

Catalog No	YP-mAb-11519
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LRRC59 PRO1855
Immunogen	Synthesized peptide derived from human LRC59 AA range: 110-160
Specificity	This antibody detects endogenous levels of LRC59 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Leucine-rich repeat-containing protein 59;N-terminally processed;p34;PRO1855;LRC59;LRRC59 PRO1855
Observed Band	34kD
Calculated Molecular Weight	34kD
Cell Pathway	Microsome membrane ; Single-pass type II membrane protein . Endoplasmic reticulum membrane; Single-pass type II membrane protein. Nucleus envelope. Localization in the nuclear envelope depends upon the nuclear import machinery, including KPNB1.
Tissue Specificity	Widely expressed.
Function	similarity:Contains 5 LRR (leucine-rich) repeats.,subunit:Interacts with SGOL1.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using LRC59 mouse mAb