



TRML4 mouse mAb

Catalog No	YP-mAb-11520
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TREML4 TLT4 UNQ9425/PRO34675
Immunogen	Synthesized peptide derived from human TRML4 AA range: 140-190
Specificity	This antibody detects endogenous levels of TRML4 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TLT 4;TLT4;Trem like transcript 4 protein;TREML4;TRML4;TREML4 TLT4 UNQ9425/PRO34675
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	
Function	similarity:Contains 1 Ig-like V-type (immunoglobulin-like) domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

Western Blot analysis of various cells using TRML4 mouse mAb

