




GNT2C mouse mAb

Catalog No	YP-mAb-11537
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GCNT2 GCNT5 II NACGT1
Immunogen	Synthesized peptide derived from human GNT2C AA range: 251-301
Specificity	This antibody detects endogenous levels of GNT2C at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CCAT;EC:2.4.1.150;GCNT2C;GCNT5;I-branching enzyme;GNT2C;GCNT2 GCNT5 II NACGT1
Observed Band	
Calculated Molecular Weight	44kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using GNT2C mouse mAb