




CA021 mouse mAb

Catalog No	YP-mAb-11582
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	C1orf21 PIG13
Immunogen	Synthesized peptide derived from human CA021 AA range: 24-74
Specificity	This antibody detects endogenous levels of CA021 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CA021;C1orf21 PIG13
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Expressed in spleen, prostate, testis and uterus.
Function	tissue specificity:Expressed in spleen, prostate, testis and uterus.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using CA021 mouse mAb