




SPT22 mouse mAb

Catalog No	YP-mAb-11613
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SPATA22
Immunogen	Synthesized peptide derived from human SPT22 AA range: 196-246
Specificity	This antibody detects endogenous levels of SPT22 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NYD SP20;SPATA 22;Spermatogenesis-associated protein 22;Testis development;SPT22;SPATA22
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Chromosome . Localizes on meiotic chromosome axes. Accumulates on resected DNA. Localization is dependent on MEIOB. .
Tissue Specificity	Highly expressed in adult testis.
Function	tissue specificity:Highly expressed in adult testis.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using SPT22 mouse mAb