



RPIA mouse mAb

Catalog No	YP-mAb-11617
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RPIA RPI
Immunogen	Synthesized peptide derived from human RPIA AA range: 256-306
Specificity	This antibody detects endogenous levels of RPIA at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RPI;Ribose-5-phosphate isomerase;ribose 5 phosphate isomerase A; Phosphoriboisomerase;EC:5.3.1.6;RPIA;RPIA RPI
Observed Band	
Calculated Molecular Weight	
Cell Pathway	cytosol,integral component of membrane,intracellular membrane-bounded organelle,
Tissue Specificity	
Function	catalytic activity:D-ribose 5-phosphate = D-ribulose 5-phosphate., disease:Defects in RPIA are the cause of ribose 5-phosphate isomerase deficiency [MIM:608611]. A patient has been described with a deficiency of ribose 5-phosphate isomerase who presented with leukoencephalopathy and peripheral neuropathy. Proton magnetic resonance spectroscopy of the brain revealed a highly elevated level of the polyols ribitol and D-arabitol, which were subsequently also found in high concentrations in body fluids. Deficient activity of RPIA, one of the pentose phosphate pathway enzymes, has been demonstrated in fibroblasts.,pathway:Carbohydrate degradation; pentose phosphate pathway; D-ribose 5-phosphate from D-ribulose 5-phosphate (non-oxidative stage): step 1/1.,similarity:Belongs to the ribose 5-phosphate



isomerase family.,

Background

The protein encoded by this gene is an enzyme, which catalyzes the reversible conversion between ribose-5-phosphate and ribulose-5-phosphate in the pentose-phosphate pathway. This gene is highly conserved in most organisms. The enzyme plays an essential role in the carbohydrate metabolism. Mutations in this gene cause ribose 5-phosphate isomerase deficiency. A pseudogene is found on chromosome 18. [provided by RefSeq, Mar 2010],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using RPIA mouse mAb