




DET1 mouse mAb

Catalog No	YP-mAb-11619
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DET1
Immunogen	Synthesized peptide derived from human DET1 AA range: 48-98
Specificity	This antibody detects endogenous levels of DET1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DET1
Observed Band	
Calculated Molecular Weight	61kD
Cell Pathway	Nucleus .
Tissue Specificity	Highly expressed in the ovary, some lymphoid organs and resting leukocytes.
Function	function:Component of the E3 ubiquitin ligase DCX DET1-COP1 complex, which is required for ubiquitination and subsequent degradation of target proteins. The complex is involved in JUN ubiquitination and degradation.,similarity:Belongs to the DET1 family.,subunit:Component of the DCX DET1-COP1 ubiquitin ligase complex at least composed of RBX1, DET1, DDB1, CUL4A and COP1.,tissue specificity:Highly expressed in the ovary, some lymphoid organs and resting leukocytes.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using DET1 mouse mAb