




WWC3 mouse mAb

Catalog No	YP-mAb-11621
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	WWC3 KIAA1280
Immunogen	Synthesized peptide derived from human WWC3 AA range: 363-413
Specificity	This antibody detects endogenous levels of WWC3 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	WWC3;WWC3 KIAA1280
Observed Band	
Calculated Molecular Weight	
Cell Pathway	cytosol,
Tissue Specificity	
Function	similarity:Belongs to the WWC family.,similarity:Contains 1 C2 domain.,
Background	This gene encodes a member of the WWC family of proteins, which also includes the WWC1 (KIBRA) gene product and the WWC2 gene product. The protein encoded by this gene includes a C2 domain, which is known to mediate homodimerization in the related WWC1 gene product. [provided by RefSeq, Sep 2011],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using WWC3 mouse mAb