




CCD70 mouse mAb

Catalog No	YP-mAb-11622
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CCDC70
Immunogen	Synthesized peptide derived from human CCD70 AA range: 149-199
Specificity	This antibody detects endogenous levels of CCD70 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CCDC70;CCD70
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Mitochondrion .
Tissue Specificity	
Function	caution:It is uncertain whether Met-1 or Met-12 is the initiator.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using CCD70 mouse mAb