




GSG1 mouse mAb

Catalog No	YP-mAb-11628
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GSG1 UNQ709/PRO1360
Immunogen	Synthesized peptide derived from human GSG1 AA range: 24-74
Specificity	This antibody detects endogenous levels of GSG1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	UNQ709/PRO1360;Germ cell-specific gene 1 protein;germ cell associated 1; GSG1;GSG1 UNQ709/PRO1360
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein . Colocalizes with PAPOLB in the endoplasmic reticulum. .
Tissue Specificity	
Function	function:May cause the redistribution of PAPOLB from the cytosol to the endoplasmic reticulum.,similarity:Belongs to the GSG1 family.,subcellular location:Colocalizes with PAPOLB in the endoplasmic reticulum., subunit:Interacts with PAPOLB.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using GSG1 mouse mAb