




ZBBX mouse mAb

Catalog No	YP-mAb-11638
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZBBX
Immunogen	Synthesized peptide derived from human ZBBX AA range: 295-345
Specificity	This antibody detects endogenous levels of ZBBX at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Zinc finger B box domain containing protein 1;ZBBX
Observed Band	
Calculated Molecular Weight	
Cell Pathway	intracellular,
Tissue Specificity	
Function	similarity:Contains 1 B box-type zinc finger.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using ZBBX mouse mAb